




MedGyn Blue Titanium

Classic Tischler bite with durable coating.
Bite Size: 7mm x 3mm

030176 Tischler Biopsy Punch Blue Titanium
Heavy-Grip Handle 10 (254mm)

030181 Tischler Biopsy Punch Blue Titanium 10
(254mm)



Baby Tischler Biopsy Punch

Smaller bite size designed for less trauma with durable coating. Bite Size: 4.2mm x 2.3mm

030177 Baby Tischler Biopsy Punch Blue Titanium
Heavy-Grip Handle 10 (254mm)

030129 Baby Tischler Biopsy Punch Blue Titanium
8½ (216mm)

030180 Baby Tischler Biopsy Punch Blue Titanium
10 (254mm)



Baby Tischler Biopsy Punch (Forceps Style)

Smaller bite size designed for less trauma with durable coating. Bite Size: 4.2mm x 2.3mm

030179 Baby Tischler Biopsy Punch Blue Titanium
Gold Handle (Forceps Style) 10" (254mm)



Kevorkian Biopsy Punch

Popular design with strong, narrow jaw with durable coating and a heavy-grip handle. Bite Size: 8mm x 3mm

030178 Kevorkian Biopsy Punch Blue Titanium
Heavy-Grip Handle 10 (254mm)



Oval Tischler Biopsy Punch

Oval Tischler bite with durable coating and a heavy-grip handle.

030175 Oval Tischler Biopsy Punch Blue Titanium
Heavy-Grip Handle 10 (254mm)



MedGyn Titanium

Tischler Biopsy Punch

Classic Tischler bite, with durable coating.

Bite Size: 7mm x 3mm

030172 Tischler Biopsy Punch Titanium 10 (254mm) Heavy-Grip Handle

030120 Tischler Biopsy Punch Titanium 10 (254mm)

030124 Tischler Biopsy Punch Titanium 10 (254mm)-Angled Up

030125 Tischler Biopsy Punch Titanium 10 (254mm)-Angled Down

030162 Tischler Biopsy Punch Titanium 10 (254mm)-Angled Left

030163 Tischler Biopsy Punch Titanium 10 (254mm)-Angled Right

030118 Tischler Biopsy Punch Titanium 8½ (216mm)

030122 Tischler Biopsy Punch Titanium 8½ (216mm)-Angled Up

030123 Tischler Biopsy Punch Titanium 8½ (216mm)-Angled Down

030138 Tischler Biopsy Punch Titanium 8½ (216mm)-Angled Right

030139 Tischler Biopsy Punch Titanium 8½ (216mm)-Angled Left



Baby Tischler Biopsy Punch (Forceps Style)

Classic Baby Tischler bite with durable coating. Bite Size:

4.2mm x 2.3mm

030170 Baby Tischler Biopsy Punch Titanium Gold Handle (Forceps Style) 10" (254mm)



Oval Tischler Biopsy Punch

030171 Oval Tischler Biopsy Punch Titanium Heavy-Grip Handle 10 (254mm)



Kevorkian Biopsy Punch

Popular design with strong, narrow jaw with durable coating.

Bite Size: 8mm x 3mm.

030174 Kevorkian Biopsy Punch Titanium Heavy-Grip Handle 10" (254mm)

030127 Kevorkian Biopsy Punch Titanium 10 (254mm)

030126 Kevorkian Biopsy Punch Titanium 8½ (216mm)



Baby Tischler Biopsy Punch

Smaller bite size designed for less trauma with durable coating. Bite Size: 4.2mm x 2.3mm

- 030173 Baby Tischler Biopsy Punch Titanium 10 (254mm)
Heavy-Grip Handle
- 030121 Baby Tischler Biopsy Punch Titanium 10 (254mm)
- 030134 Baby Tischler Biopsy Punch Titanium 10 (254mm)
- Angled Up
- 030135 Baby Tischler Biopsy Punch Titanium 10 (254mm)
- Angled Down
- 030161 Baby Tischler Biopsy Punch Titanium 10 (254mm)
- Angled Left
- 030160 Baby Tischler Biopsy Punch Titanium 10 (254mm)
- Angled Right
- 030119 Baby Tischler Biopsy Punch Titanium 8½ (216mm)
- 030132 Baby Tischler Biopsy Punch Titanium 8½ (216mm)
- Angled Up
- 030133 Baby Tischler Biopsy Punch Titanium 8½ (216mm)
- Angled Down
- 030159 Baby Tischler Biopsy Punch Titanium 8½ (216mm)
- Left
- 030158 Baby Tischler Biopsy Punch Titanium 8½ (216mm)
- Right



MedGyn Stainless Steel

Tischler Biopsy Punch

Classic Tischler bite designed for drawing in maximum tissue near Os with minimal trauma.

Bite Size 7mm x 3mm

- 030101-10 Tischler Biopsy Punch 10 (254mm)
- 029801 Tischler Biopsy Punch 10 (254mm) Angled Up
- 029802 Tischler Biopsy Punch 10 (254mm) Angled Down
- 029803 Tischler Biopsy Punch 10 (254mm) Angled Left
- 029804 Tischler Biopsy Punch 10 (254mm) Angled Right
- 030101 Tischler Biopsy Punch 8½ (216mm)
- 030113 Tischler Biopsy Punch 8½ (216mm) Angled Up
- 030114 Tischler Biopsy Punch 8½ (216mm) Angled Down
- 030136 Tischler Biopsy Punch 8½ (216mm) Angled Left
- 030137 Tischler Biopsy Punch 8½ (216mm) Angled Right



Baby Tischler Biopsy Punch

Smaller bite size designed for less trauma.

Bite Size: 4.2mm x 2.3mm

030102-10 Baby Tischler Biopsy Punch 10 (254mm)

029811 Baby Tischler Biopsy Punch 10 (254mm) Angled Up

029812 Baby Tischler Biopsy Punch 10 (254mm) Angled Down

029813 Baby Tischler Biopsy Punch 10 (254mm) Angled Left

029814 Baby Tischler Biopsy Punch 10 (254mm) Angled Right


030102 Baby Tischler Biopsy Punch 8½ (216mm)

030140 Baby Tischler Biopsy Punch 8½ (216mm) Angled Up

030141 Baby Tischler Biopsy Punch 8½ (216mm) Angled Down

030142 Baby Tischler Biopsy Punch 8½ (216mm) Angled Left

030143 Baby Tischler Biopsy Punch 8½ (216mm) Angled Right



Kevorkian Biopsy Punch

(Kevorkian-Younger Biopsy Forceps displayed)

030103-10 Kevorkian Biopsy Punch 10 (254mm)

030103 Kevorkian Biopsy Punch 8½ (216mm)

030169 Kevorkian-Younger Biopsy Forceps without Spring
8" (203mm)



Tischler-Kevorkian Cervical Biopsy Punch

030116 Tischler-Kevorkian Biopsy Punch 10
(254mm)

030115 Tischler-Kevorkian Biopsy Punch 8½
(216mm)



Burke Biopsy Punch


Open-handle design for optimal control with small cutting jaws and sharp teeth.

030109 Burke Biopsy Punch 8½



Eppendorfer Biopsy Forceps

030110 Eppendorfer Biopsy Forceps 8" (203mm)



Van Doren Biopsy Forceps

030155 Van Doren Biopsy Forceps with Ratchet 10" (254mm)

030166 Van Doren Biopsy Forceps without Ratchet 10" (254mm)



Schubert Biopsy Forceps

030157 Schubert Biopsy Forceps 10" (254mm) Curved



Wittner Biopsy Punch

030117 Wittner Biopsy Punch 8½ (216mm)



Rotating Biopsy Punch Options

030778 Vulvar Keyes Biopsy Punch Set (Handle + 6 Tips)

030104 Rotating Tischler Biopsy Punch

030145 Rotating Tischler Biopsy Punch - Angled Down

030146 Rotating Tischler Biopsy Punch - Angled Up

030105 Rotating Baby Tischler Biopsy Punch

030105-DOWN Rotating Baby Tischler Biopsy Punch - Angled Down

030105-UP Rotating Baby Tischler Biopsy Punch - Angled Up

030106 Rotating Kevorkian Biopsy Punch

030107 Rotating Burke Biopsy Punch

030111-BB Rotating Shaft - Baby Tischler Tip Only

030111-BRK Rotating Shaft - Burke Tip Only

030111-KEV Rotating Shaft - Kevorkian Tip Only

030111-TIS Rotating Shaft - Tischler Tip Only

030147 Rotating Shaft - Tischler Tip - Up

030148 Rotating Shaft - Tischler Tip - Down

030108 Rotating Biopsy Punch Handle Only

030100 Rotating Shaft (with any tip)

