

**Vortex Surgical, USA**

Supplier of high quality disposable forceps, laser probes, intravitreal procedure kits, cannulae and lenses for retinal surgery.

**Myopic ACTU8 Forceps****NOW AVAILABLE!**

Perfect for high myopic eyes with a 5mm increased shaft length

**All the same benefits of the ACTU8 regular line of forceps including:**

- 360 degree handle
- Highly responsive actuation mechanism
- Texturized finish for enhanced grip strength and glare free finish

**Myopic ACTU8 ILM Forceps**

- High visibility design
- Designed to grab and hold the ILM
- Minimises membrane stripping
- Ideal for ILM and ERM
- Sterile box of 5

**Myopic ACTU8 Tewari T2 Forceps**

- Tip designed to purchase ILM with buffering to retinal surface.
- High quality design
- Designed by Dr. Asheesh Tewari for use on ILM, ERM and tissue engagement of all types
- Sterile box of 5

- VS0744.23-35 - 23ga Myopic ILM Forceps - Box 5 sterile packed - \$957.00 + GST
- VS0744.25-35 - 25ga Myopic ILM Forceps - Box 5 sterile packed - \$957.00 + GST
- VS0742.23-35 - 23ga Myopic Tewari T2 Forceps - Box 5 sterile packed - \$957.00 + GST
- VS0742.25-35 - 25ga Myopic Tewari T2 Forceps - Box 5 sterile packed - \$957.00 + GST