


Distributed in New Zealand by:




Cooper Medical Limited
57 Mahunga Drive
Mangere Bridge
Auckland 2022

PMS



BASIC



GENERAL INSTRUMENTS 16

SPECULAE & RETRACTORS 17

BLADE BREAKERS, HANDLES & BLADES 18

CALIPERS & FIXATION RINGS 19

CANNULAE & ACCESSORIES 20

IOL & PHACO INSTRUMENTS 21

SCISSORS 22

FORCEPS +
TYPES OF FORCEPS TIPS 23-24

NEEDLE HOLDERS 25

MISCELLANEOUS PRODUCTS 26


DISPOSABLE PRODUCTS 27

BASIC EYE SET 28






GENERAL INSTRUMENTS




**Halstead
Mosquito Forceps**
12.5 cm

- *E30-902 straight
- *E30-903 curved





**Dieffenbach
Serrefine, 3.5 cm**

- E30-934 straight
- E30-935 curved



**Baby Jones
Towel Clamp, 5.5 cm**

- E30-965




**Knapp
Strabismus Scissors**
short blades, blunt, 12.0 cm

- E25-556 straight
- E25-557 curved



**Standard
Eye Scissors, pointed**

- E25-510 straight, 9.0 cm
- E25-512 straight, 10.5 cm
- E25-514 straight, 11.5 cm
- E25-511 curved, 9.0 cm
- E25-513 curved, 10.5 cm
- E25-515 curved, 11.5 cm




*also available in titanium, please add "T" when ordering



SPECULAE & RETRACTORS

BASIC EYE INSTRUMENTS




**Barraquer
Wire Speculum**

E40-011 adult size

E40-012 child size

E40-013 infant size

E40-014 premature size




**Kratz-Barraquer
Wire Speculum**

E40-031 adult size

E40-032 child size

E40-033 adult size,
heavy wire model



Desmarres

Lid Retractor, 14.0 cm

*E40-110 #0 11 mm blade

*E40-111 #1 13 mm blade

E40-112 #2 15 mm blade

Conway
pediatric size, 8 mm blade

*E40-108

Castroviejo

Universal Eye Speculum
fenestrated blades, 9.0 cm

E40-078 medium, 16 mm

E40-079 small, 14 mm



Fishkind
Phaco Eye Speculum
open blades, 8.5 cm

E40-080




for use with Phaco
or I/A tips

*also available on round titanium handle, please add "T" when ordering



BLADE BREAKER, HANDLES & BLADES

BASIC EYE INSTRUMENTS



Razor Blade
easy to break,
stainless steel,
non-sterile,
(box of 10)
E35-001



Troutman-Chris
Blade Breaker & Holder
5.5 mm dia. x 9.5 cm long handle
straight jaws with concave/
convex inside surfaces
***E35-010**



Chuck Handle
(for miniature edged blades)
***E10-003K #3K, 10.0 cm**
***E10-003KM #3KM, 7.5 cm**



Blade Collet
***E10-003BC**



Miniature Edged Blades
surgical stainless steel,
single use, sterile (box of 10)
E10-064 #64
E10-065 #65
E10-067 #67
E10-069 #69

**for precise
microsurgical incisions**

*also available in titanium, please add "T" when ordering

CALIPER & FIXATION RINGS

PMS




Castroviejo

Eye Caliper
0-20 mm scale
on either side in
1 mm increments

*E40-554 straight

E40-555 curved



Fine-Thornton

Phaco Fixation Ring

3/4 open swivel ring,
with multiple teeth, 13.5 cm

*E10-513 13 mm ID
for clear cornea incision

E10-516 16 mm ID
for scleral tunnel incision



Schocket

Scleral Depressor

with pocket clip, 13.5 cm

*E40-613



Flieringa

Scleral Fixation Ring

E40-628 Set of 8 sizes

E40-614 14 mm dia.

E40-615 15 mm dia.

E40-616 16 mm dia.

E40-617 17 mm dia.

E40-618 18 mm dia.

E40-619 19 mm dia.

E40-620 20 mm dia.

E40-622 22 mm dia.



Maloney

Astigmatic Keratometer

E10-690T

for qualitative evaluation
of astigmatism




*also available in titanium, please add "T" when ordering

BASIC EYE INSTRUMENTS




CANNULAE & ACCESSORIES




**Bishop-Harmon
AC Irrigator System**
E20-001 complete
E20-002 silicone bulb
E20-003 LL adaptor
E20-004 19GA cannula



**Lewicky
Anterior Chamber Maintainer**
 20GA/0.9 mm x 3 mm long screw threaded tip
 attached to 180 mm silicone tubing with LL cone
E20-020ST




**Rycroft
Air Injection Cannula**
 20 mm tube, angled 5 mm
E20-030 30GA/0.3 mm
E20-031 27GA/0.4 mm
E20-032 25GA/0.5 mm



**Simcoe
Double I and A Cannula**
 23/23GA x 15 mm long tube, 0.3 mm aspirating port
 attached to 150 mm silicone tubing with LL adaptor
E20-300 original model
E20-301 reverse model

All Glass Syringe **E20-801** 1 ml
E20-802 2 ml
E20-805 5 ml

ideal for smooth control of
irrigation or maintenance
of the anterior chamber




**McIntyre
Aspirating Handpiece**
 (autoclavable) Teflon, 7.5 cm
 male/male Luer connectors
E20-881

Luer-Lock cone

IOL & PHACO INSTRUMENTS



Jaffe-Bechert
Lens Nucleus Rotator
angled 9 mm shaft,
blunt forked tip,
11.5 cm
*E15-423



Sinsky II
IOL Manipulating Hook
0.2 mm dia. blunt tip,
12.0 cm
*E15-552 straight
*E15-553 angled



Lester
Lens Pusher
delicate "hour glass" style tip,
12.0 cm
*E15-500 straight
*E15-501 angled



Kuglen
Lens Manipulator
and Iris Hook
"push-pull" model,
angled 10 mm shaft,
ball-like tip, 12.0 cm



*E15-565

Nagahara
Phaco Chopper
angled 10 mm shaft,
"wedge-shaped" blade,
11.5 cm
*E15-656 60° cutting edge
*E15-659 90° cutting edge



introduced in 1975 still
as popular as ever



also available on coloured aluminium handles, see page 95, 107 and 131

*also available on titanium handle, see page 34-36


BASIC EYE INSTRUMENTS




SCISSORS




**Barraquer
Iris Scissors**
7 mm blades,
blunt tips, 5.5 cm
E25-200



Micro Iris Scissors
curved, small blades,
10.0 cm
*E25-251 sharp tips
*E25-253 blunt tips




**Westcott
Stitch Scissors**
curved, large blades,
sharp tips, 11.5 cm
*E25-315



**the working horse
of cataract surgery**



**Westcott
Tenotomy Scissors**
curved, large blades,
blunt tips, 11.5 cm
E25-334T right (shown)
E25-334 right
E25-335 left



**Vannas
Capsulotomy Scissors**
5 mm blades, 8.2 cm
E25-400 straight, sharp
E25-401 curved, sharp
E25-402 straight, blunt
E25-403 curved, blunt


all time favourite

*also available on titanium handle, please add "T" when ordering




FORCEPS

TYPES OF FORCEPS TIPS




1x2 teeth 90°
(Bonn style)



1x2 teeth 45°
(Castroviejo style)



serrated jaws



toothed jaws



criss-cross serrated



Barraquer-Colibri
Corneal Utility Forceps
1x2 teeth 90°, 0.12 mm,
without tying platform, 7.5 cm

*E30-051



Harms-Colibri
Corneal Suture Forceps
1x2 teeth 45°, 0.12 mm,
tying platform, 7.5 cm

*E30-081



Dressing Forceps
serrated jaws, 10.0 cm
E30-110 0.5 mm wide
E30-114 1.0 mm wide



Tissue Forceps
toothed jaws, 10.0 cm
E30-120 0.5 mm teeth
E30-124 1.0 mm teeth



Bishop-Harmon
Dressing & Tissue Forceps
E30-100 criss-cross serrated
E30-103 0.3 mm teeth
E30-106 0.6 mm teeth

*also available in titanium, see page 44








FORCEPS

these forceps are ideal...

...for 11-0 to 8-0 sutures

still the most popular



Bonn
Suturing Forceps
1x2 teeth 90°,
0.12 mm, delicate
4 mm platform
*E30-410 short model,
7.2 cm

*E30-430 long model, 9.5 cm



Harms
Tying Forceps
6 mm platform, 10.0 cm

*E30-422 straight

*E30-423 curved



McPherson
Tying Forceps
5 mm platform, 10.5 cm

*E30-460 straight

*E30-461 angled



Castroviejo (standard handle)
Suturing Forceps
1x2 teeth 45°, fine
tying platform, 10.0 cm

*E30-431 0.12 mm

E30-432 0.30 mm

*E30-433 0.50 mm



Utrata
Capsulorhexis Forceps
extra thin 11 mm shanks,
triangular, sharp tips, 10.5 cm

*E30-581S flat handle

*E30-583S knurled handle


*also available in titanium, see page 41,43,45

NEEDLE HOLDERS

PMS




straight jaws
E35-130
without lock
*E35-132
with lock



**Castroviejo
Needle Holders**
14.0 cm

delicate jaws, 10 mm x 1.0 mm
ideal for 7-0 to 4-0 sutures

curved jaws
E35-131
without lock
*E35-133
with lock



straight jaws
E35-230
without lock
E35-232
with lock



**Barraquer
Needle Holders**
13.5 cm

delicate jaws, 10 mm x 0.75 mm
ideal for 9-0 to 6-0 sutures

curved jaws
*E35-231
without lock
*E35-233
with lock



*also available in titanium, please add "T" when ordering

BASIC EYE INSTRUMENTS


ECED0308F



MISCELLANEOUS PRODUCTS



Protector
for Lens Hooks with
E15-002 straight shaft
E15-003 angled shaft



Metal Ruler
flexible, graduated in mm & inches
E40-559

RS-SPRAY
Rust Prevention Spray (sterilisable)
for stainless steel instruments,
spray box of 200 ml
No. 2820

Ideal to prevent rusting of stainless steel surgical instruments during and after sterilisation. Applicable for all methods of sterilisation. Autoclave max. 134°C and hot air sterilizers max. 180°C.

Type	Id. (mm)	Wd. (mm)	Od. (mm)
Product No.			
SIL-00301	0,3	0,1	0,5
SIL-00302	0,3	0,2	0,7
SIL-00502	0,5	0,2	0,9
SIL-00503	0,5	0,3	1,1
SIL-00504	0,5	0,4	1,3
SIL-00505	0,5	0,5	1,5
SIL-00701	0,7	0,1	0,9
SIL-00702	0,7	0,2	1,1
SIL-00703	0,7	0,3	1,3
SIL-00705	0,7	0,5	1,7
SIL-00804	0,8	0,4	1,6
SIL-01002	1,0	0,2	1,4
SIL-01004	1,0	0,4	1,8
SIL-01005	1,0	0,5	2,0
SIL-01010	1,0	1,0	3,0
SIL-01204	1,2	0,4	2,0
SIL-01502	1,5	0,2	1,9
SIL-01503	1,5	0,3	2,1
SIL-01504	1,5	0,4	2,3
SIL-01505	1,5	0,5	2,5
SIL-01507	1,5	0,7	2,9
SIL-01510	1,5	1,0	3,5
SIL-02005	2,0	0,5	3,0
SIL-02007	2,0	0,7	2,9
SIL-02010	2,0	1,0	4,0
SIL-02505	2,5	0,5	3,5
SIL-02510	2,5	1,0	4,5

Silicone Tubing (VMQ Quality)
autoclavable up to max. 140°C,
available in coils of 100 cm

SIL-00302 for E40-900 and E40-910

SIL-00804 for E20-015 to E20-019

SIL-01004 E20-300 to E20-319 and Lieberman style aspirating speculae (side connectors)

SIL-01204 for E20-321 and E20-323

SIL-01504 for E20-751 to E20-753 and Lieberman aspirating speculae (main connectors)

SIL-02005 for E20-741, E20-754 to E20-757

Further sizes are available by request and in coils of 25 metres.



DISPOSABLE PRODUCTS



Sugi®
Cellulose Eye Spears
white plastic handle,
single use, sterile
30601 box of 200
(20 pouches of 10 pcs)

highly absorbent sponge for removal of blood,
fluids and cellular debris within the operating field

Sugi®
Cellulose Eye Points
single use, sterile
31602 box of 200 (20 pouches of 10 pcs)
31603 box of 500 (bulk, non-sterile)

Pro-ophta®
Ocular Sticks
5 mm dia. x 66 mm,
single use, sterile
16516 box of 400 (20 pouches of 20 pcs)
16515 box of 500 (bulk, non-sterile)

Insight®
Universal Eye Shield
single, sterile blister packed (box of 100)
ES/100/C

dual fitting to either eye
with ventilation holes for vision

Universal Eye Shield
single, sterile blister packed (box of 100)
ES/100/W



crystal clear K Resine material

Latex free

white polypropylene material

BASIC EYE INSTRUMENTS



E40-086



E30-232



E30-532



E40-771



E40-777

- E20-001 AC Irrigator System, page 20
- E20-031 Air Injection Cannula, page 20
- E30-965 Schaedel Towel Clamp, page 16
- E30-902 Halstead Artery Forceps, page 16
- E30-934 Dieffenbach Serrefine, page 16
- E40-011 Barraquer Eye Speculum, page 17
- E40-086 Lancaster Universal Eye Speculum**
- E10-014 Sapphire Knife 45° blade, page 63
- E35-001 Razor Blade (box of 10), page 18
- E35-010 Blade Breaker & Holder, page 18
- E25-200 DeWecker Iris Scissors, page 22
- E25-315 Westcott Stitch Scissors, page 22
- E25-334 Westcott Tenotomy Scissors, page 22
- E25-401 Vannas Capsulotomy Scissors, page 22
- E25-512 Standard Eye Scissors, page 16
- E25-556 Knapp Strabismus Scissors, page 16
- E30-051 Barraquer Colibri Forceps, page 23
- E30-106 Bishop-Harmon Tissue Forceps, page 23
- E30-110 Dressing Forceps, straight, page 23
- E30-232 Lester Conjunctival Fixation Forceps**
- E30-431 Castroviejo Suture Forceps, page 24
- E30-432 Castroviejo Suture Forceps, page 24
- E30-460 McPherson Tying Forceps, page 24
- E30-463 Kelman-McPherson Forceps, page 124
- E30-532 Bipolar Forceps "Jewelers" Style**
- E35-231 Barraquer Needle Holder, page 25
- E35-132 Castroviejo Needle Holder, page 25
- E40-111 Desmarres Lid Retractor, page 17
- E40-559 Stainless Steel Ruler, page 26
- E40-554 Castroviejo Eye Caliper, page 19
- E40-613 Schocket Scleral Depressor, page 19
- E40-771 Graefe Strabismus Hook**
- E40-777 Gass Retinal Detachment Hook**
- Schioetz Tonometer, Wash Trays and Sterilising Trays (see section)