

Coaxial Disposable Irrigation Aspiration System 20G, 21G, 22G




Product Features:

- Coaxial System Disposable 20G & 22G, 21G optional
- Lightweight irrigation / aspiration Handpiece
- Prolonged sleeve (plus 1.5mm)
- all G: 45° tip smooth abrasive (rough) or polished
- 20G, 21G: 65° tip smooth abrasive (rough) or polished
- Compatible with your existing Phaco or I/A System
- Swiss Design / Swiss Made


Details:




prolonged tip



45° tip



65° tip



./.

Distributed in New Zealand by:



Cooper Medical Limited
57 Mahunga Drive
Mangere Bridge
Auckland 2022

About Bürki
Bürki plans, develops and produces
microsurgical instruments in SWISS QUALITY
EN ISO 13485

Bürki inno med AG
Industriestrasse 67
9443 Widnau
Switzerland

Tel +41 (0)71 744 88 66
info@buerki.eu
www.buerki.eu

CX2045-P

Coaxial System Disposable 20G
Asp. port 0.35mm, angled 45°, POLISHED

CX2065-P

Coaxial System Disposable 20G
Asp. port 0.35mm, angled 65°, POLISHED

CX2145-P

Coaxial System Disposable 21G
Asp. port 0.35mm, angled 45°, POLISHED

CX2165-P

Coaxial System Disposable 21G
Asp. port 0.35mm, angled 65°, POLISHED

CX2245-P

Coaxial Sytem Disposable 22G
Asp. port 0.35mm, angled 45°, POLISHED

CX2045-R

Coaxial System Disposable 20G
Asp. Port 0.35mm, angled 45°, SMOOTH
ABRASIVE (ROUGH)

CX2065-R

Coaxial System Disposable 20G
Asp. port 0.35mm, angled 65°, SMOOTH
ABRASIVE (ROUGH)

CX2145-R

Coaxial System Disposable 21G
Asp. port 0.35mm, angled 45°, SMOOTH
ABRASIVE (ROUGH)

CX2165-R

Coaxial System Disposable 21G
Asp. port 0.35mm, angled 65°, SMOOTH
ABRASIVE (ROUGH)

CX2245-R

Coaxial Sytem Disposable 22G
Asp. port 0.35mm, angled 45°, SMOOTH
ABRASIVE (ROUGH)

Version 05/2022

About Bürki
Bürki plans, develops and produces
microsurgical instruments in SWISS QUALITY
EN ISO 13485

Bürki inno med AG
Industriestrasse 67
9443 Widnau
Switzerland

Tel +41 (0)71 744 88 66
info@buerki.eu
www.buerki.eu