

Bimanual System Disposable Irrigation Aspiration 22G, 23G



Distributed in New Zealand by:



Cooper Medical Limited
57 Mahunga Drive
Mangere Bridge
Auckland 2022

22G



BSD 540

Bimanual System Disposable **22G** consisting of:

1 Irrigation Handpiece, with curved cannula and dual side ports
ø 0.45mm

1 Aspiration Handpiece with curved cannula and asp. port ø 0.35 mm,
smooth abrasive (rough) tip



BSD 541

Bimanual System Disposable **22G** consisting of:

1 Irrigation Handpiece, with curved cannula and dual side ports
ø 0.45 mm

1 Aspiration Handpiece with curved cannula and asp. port ø 0.35 mm,
polished tip



BSD 542

Bimanual System Disposable **22G** consisting of:

1 Irrigation Handpiece, with curved cannula and open end flushing 45°
polished

1 Aspiration Handpiece with curved cannula and asp. port ø 0.35 mm,
smooth abrasive (rough) tip

23G



BSD 522

Bimanual System Disposable **23G** consisting of:

1 Irrigation Handpiece with curved cannula and open end flushing 45°,
polished

1 Aspiration Handpiece with curved cannula and asp. port ø 0.30 mm,
smooth abrasive (rough) tip

BSD 520 and BSD 521 soon available



BSD 520

Bimanual System Disposable **23G** consisting of:

1 Irrigation Handpiece, with curved cannula, dual side ports ø 0.35 mm

1 Aspiration Handpiece, with curved cannula, Asp. Port ø 0.30 mm
smooth abrasive (rough) tip



BSD 521

Bimanual System Disposable **23G** consisting of:

1 Irrigation Handpiece with curved cannula, dual side ports ø 0.35 mm

1 Aspiration Handpiece with curved cannula, Asp. Port ø 0.30 polished
tip

About Bürki
Bürki plans, develops, produces and certifies
microsurgical instruments in SWISS QUALITY
DIN EN ISO 13485

Bürki inno med AG
Industriestrasse 67
9443 Widnau
Switzerland

Tel +41 (0)71 744 88 66
Fax +41 (0)71 744 88 67
info@buerki.eu
www.buerki.eu

Detailed product Features:

	BSD 540	BSD 541	BSD 542	BSD 520	BSD 521	BSD 522
Single use / no risk of cross contamination	●	●	●	●	●	●
Lightweight, colour-coded handpiece	●	●	●	●	●	●
Standard male/female Luer connector	●	●	●	●	●	●
Compatible with your Phaco or I/A System	●	●	●	●	●	●
Smooth bullet-shaped tip		●			●	
Irrigation front flushing 45°			●			●
Irrigation dual side ports 0.35 mm				●	●	
Irrigation dual side ports 0.45 mm	●	●				
Aspiration port 0.30 mm, rough tip				●		●
Aspiration port 0.30 mm, smooth tip					●	
Aspiration port 0.35 mm, rough tip	●		●			
Aspiration port 0.35 mm, smooth tip		●				

Version 04/2022

About Bürki
Bürki plans, develops, produces and certifies
microsurgical instruments in SWISS QUALITY
DIN EN ISO 13485

Bürki inno med AG
Industriestrasse 67
9443 Widnau
Switzerland

Tel +41 (0)71 744 88 66
Fax +41 (0)71 744 88 67
info@buerki.eu
www.buerki.eu