

Reusable Bimanual Irrigation Aspiration System 22G, 23G



Product Features:

- Lightweight, colour coded irrigation and aspiration Handpieces
- Compatible with your existing Phaco or I/A System
- Smooth abrasive (rough) tip available for capsule polishing
- Swiss Design / Swiss Made

Distributed in New Zealand by:



Cooper Medical Limited
57 Mahunga Drive
Mangere Bridge
Auckland 2022

Aspiration

Irrigation

22G



BA 118/T
Aspiration Handpiece, polished tip,
(0.70mm), Asp. port 0.35 mm



BI 128/T Special model
Irrigation Handpiece, polished tip,
dual side ports 0.45 mm




BA 119/T Special model
Aspiration Handpiece, smooth
abrasive (rough) tip, (0.70mm)
Asp. port 0.35 mm



BI 129/T Special model
Irrigation Handpiece, polished tip,
open end Tip 45°

23G



BA 113/T
Aspiration Handpiece, polished tip,
(0.65mm) Asp. port 0.30 mm



BI 123/T
Irrigation Handpiece, round cannula,
polished tip, dual 0.35 mm side ports



BA 114/T
Aspiration Handpiece, smooth
abrasive (rough) tip,(0.65mm)
Asp. port 0.30 mm



BI 124/T
Irrigation Handpiece, polished tip
open end Tip 45° (0.65mm)

Version 04/2022

About Bürki
Bürki plans, develops, produces and certifies
microsurgical instruments in SWISS QUALITY
EN ISO 13485

Bürki inno med AG
Industriestrasse 67
9443 Widnau
Switzerland

Tel +41 (0)71 744 88 66
Fax +41 (0)71 744 88 67
info@buerki.eu
www.buerki.eu