


Bimanual System Disposable Irrigation Aspiration 21G, 20G



21G



BSD 500

Bimanual System Disposable **21G** consisting of:
1 Irrigation Handpiece with curved cannula, open end flushing 45° tip
1 Aspiration Handpiece with curved cannula, Asp. Port ø 0.35 mm
smooth abrasive (rough) tip


BSD 510

Bimanual System Disposable **21G** consisting of:
1 Irrigation Handpiece with curved cannula, dual side ports ø 0.45 mm
1 Aspiration Handpiece with curved cannula, Asp. Port ø 0.35 mm
smooth abrasive (rough) tip

BSD 515

Bimanual System Disposable **21G** consisting of:
1 Irrigation Handpiece with curved cannula, dual side ports ø 0.45 mm
1 Aspiration Handpiece with curved cannula, Asp. Port ø 0.35 mm
polished tip

20G



BSD 560

Bimanual System Disposable **20G** consisting of:
1 Irrigation Handpiece with curved cannula, dual side ports ø 0.65 mm
1 Aspiration Handpiece with curved cannula, Asp. Port ø 0.35 mm
smooth abrasive (rough) tip

BSD 565

Bimanual System Disposable **20G** consisting of:
1 Irrigation Handpiece with curved cannula, dual side ports ø 0.65 mm
1 Aspiration Handpiece with curved cannula, Asp. Port ø 0.35 mm
polished tip

Product Features:

	BSD 500	BSD 510	BSD 515	BSD 560	BSD 565
Single use / no risk of cross contamination	•	•	•	•	•
Lightweight, colour-coded Handpieces	•	•	•	•	•
Standard male/female Luer connector	•	•	•	•	•
Compatible with your Phaco or I/A System	•	•	•	•	•
Aspiration port ø 0.35 mm, roughened tip	•	•		•	
Aspiration port ø 0.35 mm, smooth tip			•		•
Irrigation front flushing 45°, inner lumen 0.6 mm	•				
Irrigation dual side ports ø 0.45 mm		•	•		
Irrigation dual side ports ø 0.65 mm				•	•
Smooth bullet-shaped tip		•	•	•	•

Version 04/2022

About Bürki
Bürki plans, develops and produces
microsurgical instruments in SWISS QUALITY
EN ISO 13485

Bürki inno med AG
Industriestrasse 67
9443 Widnau
Switzerland

Tel +41 (0)71 744 88 66
info@buerki.eu
www.buerki.eu