



010466SA-HD

010466-HD

## MedGyn AL-106 HD™ Digital Video Colposcopes

Different from the traditional optical colposcopes, the AL-106HD™ is a digital colposcope with an LED light source for ultra-bright visualization and a long life cycle. Its unique optical precision and multiple green light filters ensure the highest quality images and improved visualization of blood vessels.

- High-resolution: 1.43M px for AL-106HD, 440K for AL-106
- LED light
- 1-40x magnification
- Ultra-fast auto focus and manual focus option
- Video and image freeze
- Green light filters

- 010466-HD MedGyn AL-106 High Definition Digital Video Colposcope\*
- 010466SA-HD MedGyn AL-106 High Definition Digital Video Colposcope with Swing Arm\*
- 010010 Colposcope Monitor (Used for AL-101, AL-102, AL-401, AL-105, AL-106)\*
- 010490 Colposcope Monitor Stand (Adjustable with Rolling Base)

\*Monitor Sold Separately



010471



010469



AL-106B



AL-106 BWL

## MedGyn AL-106B and AL-106BWL Digital Video Colposcopes

The MedGyn AL-106B™ and AL-106 BWL are more affordable options to the AL-106 HD. The AL-106 BWL is a battery operated colposcope useful in remote areas with limited access to electricity.

- Battery operated and wireless
- High-resolution of 440K pixels
- LED light
- 1-54x magnification (36x for AL-106 BWL)
- Ultra-fast auto focus and manual focus option
- Video and image freeze
- Green light filters

010471 MedGyn AL-106 BWL Digital Video Colposcope, Battery Operated

010469 MedGyn AL-106B Digital Video Colposcope

\*Monitor Sold Separately