

## Reusable Bimanual Irrigation Aspiration System



### Product Features:

- Lightweight, colour coded irrigation and aspiration Handpieces
- Compatible with your existing Phaco or I/A System
- Rough tips available for capsule polishing
- Swiss Design / Swiss Made

### Most often used Bimanual I/A Systems


#### Aspiration:




**BA 100/T**  
Aspiration Handle polished tip,  
21G (0.80mm); Asp. Port  
0.30mm



**BA 102/T**  
Aspiration Handle polished tip,  
21G (0.80mm); Asp. Port  
0.40mm



**BA 104/T**  
Aspiration Handle polished tip,  
21G (0.80mm); Asp. Port  
0.50mm




**BA 105/T**  
Aspiration Handle rough tip  
21G (0.80mm); Asp. Port  
0.30mm

#### Irrigation:



**BI 120/T**  
Irrigation Handle polished tip,  
21G (0.80mm), dual side ports  
0.50mm



**BI 125/T**  
Irrigation Handle rough tip,  
21G (0.80mm), dual side ports  
0.50mm




**BI 140/T**  
Irrigation Handle polished tip,  
21G (0.80mm), open end 45°



**BI 121/T**  
Irrigation Handle polished tip,  
20G (0.90mm), dual side ports  
0.65mm

./.



**BA 108/T**  
Aspiration Handle rough tip  
20G, (0.90mm); Asp. Port  
0.35mm



**BI 126/T**  
Irrigation Handle open end tip  
45°, 20G (0.90mm), High flow



**BAO 105/T**  
Aspiration Handle oval  
cannula, 0.55mm x 0.95mm,  
mate finish, rough tip; Asp.  
Port 0.30mm



**BIO 120/T**  
Irrigation Handle oval cannula,  
0.55mm x 0.95mm, mate finish  
tip, Irrig. Port 0.55mm



**BAO 102/T**  
Aspiration Handle oval  
cannula, 0.70mm x 1.10mm,  
mate finish tip; Asp. Port  
0.35mm



**BIO 127/T**  
Irrigation Handle oval cannula,  
0.90mm x 0.75mm, mate  
finish, open end tip 30°, High  
flow