

# PRE-CLEANING



FAMOSONIC S

FAMOSONIC M

FAMOSONIC XL

Famos has various high power ultrasonic cleaning units for aqueous cleaning solutions in its product range. The benchtop models are available in three different sizes, S, M and XL. The tank of the S and M models are made of AISI 304 stainless steel, while the tank of the XL is made of AISI 316 stainless steel. All baths are equipped with a digital display including degas function.

The FamoSonic S is suitable for 1/2 DIN instrument trays such as the 50.100.063.

The FamoSonic M is suitable for 1/1 DIN instrument trays such as the 50.100.062.

The FamoSonic XL is suitable for 1/1 DIN and 1/1 ISO instrument trays such as the 50.100.062.

The ultrasonic cleaning baths need to be used in combination with an insert basket and lid.

These are available as accessories.

## BENCHTOP ULTRASONIC CLEANERS

Art.-No.	Product	Description	Packed
24.901.301	FamoSonic S	FamoSonic S benchtop ultrasonic cleaner, 13,5 liters max. capacity Internal dimensions: 325x300x150 mm LxWxH	1
24.901.302	FamoSonic M	FamoSonic M benchtop ultrasonic cleaner, 28 liters max. capacity Internal dimensions: 500x300x200 mm LxWxH	1
24.901.303	FamoSonic XL	FamoSonic XL benchtop ultrasonic cleaner, 90 liters max. capacity Internal dimensions: 600x500x300 mm LxWxH	1



Distributed in New Zealand by:

FAMOSONIC M WITH BASKET HOLDER



Cooper Medical Limited  
57 Mahunga Drive  
Mangere Bridge  
Auckland 2022

Freephone: 0800 800823 - Email: [admin@coopermedical.co.nz](mailto:admin@coopermedical.co.nz) - [www.coopermedical.co.nz](http://www.coopermedical.co.nz)

## ACCESSORIES FOR FAMOSONIC BENCHTOP ULTRASONIC CLEANERS

- Different lids available for basket holders and insert baskets
- Basket holders for 1/2 DIN, 1/1 DIN and ISO trays
- Insert baskets



In order to make sure that the ultrasonic cleaning bath works properly, it is necessary to use it in combination with a basket and lid.

For the FamoSonic S, M and XL there is an insert basket available. Alternatively, it is possible to use 1/1 DIN or ISO and 1/2 DIN trays. A tray must never be placed directly onto the bottom of the ultrasonic cleaner, because this might damage the unit. Therefore, Famos offers basket holders which the instrument trays can be placed on.



### ACCESSORIES FOR BENCHTOP ULTRASONIC CLEANERS

Art. No.	Product	Description	Packed
24.901.305	Basket holder S	Basket holder for FamoSonic S	1
24.901.306	Basket holder M / BS	Basket holder for FamoSonic M / BS	1
50.100.062	1/1 DIN tray	1/1 DIN instrument tray, 480x250x60 mm LxWxH	1
50.100.063	1/2 DIN tray	1/2 DIN instrument tray, 240x250x60 mm LxWxH	1
24.901.315	Insert basket S	Insert basket for FamoSonic S, 275x245x50 mm LxWxH	1
24.901.320	Insert basket M	Insert basket for FamoSonic M, 455x245x50 mm LxWxH	1
24.901.314	Insert basket XL	Insert basket for FamoSonic XL, 545x450x50 mm LxWxH	1
24.901.316	Lid EHK-S	Lid for FamoSonic S, only for use with insert basket	1
24.901.321	Lid EHK-M	Lid for FamoSonic M, only for use with insert basket	1
24.901.310	Lid EHK-XL	Lid for FamoSonic XL, only for use with insert basket	1
24.901.308	Lid KT-S	Lid for FamoSonic S, only for use with basket holder	1
24.901.309	Lid KT-M	Lid for FamoSonic M, only for use with basket holder	1