

## Bimanual Reusable Family Irrigation 20G, 21G, 22G, 23G

Distributed in New Zealand by:



Cooper Medical Limited  
57 Mahunga Drive  
Mangere Bridge  
Auckland 2022

**Product Features:**

- Lightweight, colour coded irrigation Handpieces
- Compatible with your existing Phaco or I/A System
- Smooth abrasive (rough) tip available for capsule polishing
- Swiss Design / Swiss Made



**Irrigation**

**20G**



**BI 126/T**  
Irrigation Handle open end tip  
45° (0.90mm), High flow

**21G**



**BI 120/T**  
Irrigation Handle polished tip  
(0.80mm), dual side ports  
0.50mm



**BI 125/T**  
Irrigation Handle smooth  
abrasive (rough) tip,  
(0.80mm), dual side ports  
0.50mm



**BI 140/T**  
Irrigation Handle polished tip  
open end 45°

**22G**



**BI 128/T special model**  
Irrigation Handle polished tip,  
dual side ports 0.45mm



**BI 129/T special model**  
Irrigation Handle polished tip,  
open end 45°

**23G**



**BI 123/T**  
Irrigation Handle polished tip  
(0.65mm), dual side ports  
0.35mm



**BI 124/T**  
Irrigation Handle polished tip  
(0.65mm), open end 45°



**BI 127/T special model (MLI)**  
Irrigation Handle polished tip  
(0.65mm)

Version 05/2022

About Bürki  
Bürki plans, develops and produces  
microsurgical instruments in SWISS QUALITY  
EN ISO 13485

Bürki inno med AG  
Industriestrasse 67  
9443 Widnau  
Switzerland

Tel +41 (0)71 744 88 66  
info@buerki.eu  
www.buerki.eu