




MASTEL SURGICAL

Where Quality Meets Precision

Interoperative markers



Distributed in New Zealand by:




Cooper Medical Limited
57 Mahunga Drive
Mangere Bridge, Auckland

Ph: 0800 800 823
sales@coopermedical.co.nz
www.coopermedical.co.nz

Handle Styles

Thornton

The traditional Thornton handle's exaggerated fenestration provides a firm grip while maintaining the traditional shape and function.



Bores

The Bores handle is a favorite of surgeons for its unique grip and fine taper. Its elegant styling also delivers a very functional nuance in the hand.

Bestselling Intraoperative Markers



- The gold standard for intraoperative Toric lens marking
- Laser etching upon an advanced coating provides detailed 5 deg marks for enhanced accuracy

Also Available:

- 25.211.0902 Gimbel Mendez 90° Thornton
- 25.211.1802 Gimbel Mendez 0/180° Thornton
- 25.211.1351 Gimbel Mendez Koch 135° Bores
- 25.211.1352 Gimbel Mendez Koch 135° Thornton
- 25.211.1001 Gimbel Mendez Jarvi Bores
- 25.211.1002 Gimbel Mendez Jarvi Thornton



- Simplifies pre-marking a true horizon by creating peripheral marks at 3 and 9 o'clock
- 3-Ray enables 3, 6 and 9 o'clock marking



- Intraoperative primary (steep) meridian marker used in conjunction with a Mendez-type Ring

Also Available: [22.102.2002](#)



- Provides accurate, fine pre-marks for astigmatic alignment that can also serve as the intraoperative toric alignment marks
- Rotating head of this plumbob handle is laser engraved in 5° increments for high levels of accuracy

Also Available:

[21.114.6025](#)

[21.114.6030](#)



- Effective in accurately marking the corneal epithelium for capsulohexis placement, or as an alcohol well in PRK surgery due to atraumatic edges and precise sizing
- Available with or without crosshairs

Also Available:

Other sizes in .25mm increments